

Appendix A – Related product codes

The product code is composed of 10 digits and explained as follows:

1	2	3	4	5	6	7	8	9	10
S	O	Obturator for I.V. catheter							
Legal manufacturer		+	Terumo Europe N.V.						
Product Type			D	M					
Indication of gauge of catheter with which the obturator is used				1	4	14 G (= 2.00 mm)			
				1	6	16 G (= 1.75 mm)			
				1	7	17 G (= 1.50 mm)			
				1	8	18 G (= 1.30 mm)			
				2	0	20 G (= 1.00 mm)			
				2	2	22 G (= 0.90 mm)			
				2	4	24 G (= 0.74 mm)			
Indication of nominal length of catheter with which the obturator is used						1	9	19 mm	
						2	5	25 mm	
						3	2	32 mm	
						4	5	45 mm	
Type of catheter with which the obturator is used								W	Winged or straight
								P	Winged and ported