

Certificate of Conformity

VWR Catalogue Number	See List
Description	Slide Boxes with Lid
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2021-08-24

Cork Lined

Catalogue Number	Colour	No. of Slides	Dimension [mm]
631-1506	Blue	25	141 × 88 × 35
631-1507	Green	25	141 × 88 × 35
631-1508	Red	25	141 × 88 × 35
631-1510	White	25	141 × 88 × 35
612-2706	Blue	50	208 × 175 × 34
612-2707	Red	50	208 × 175 × 34
612-2708	White	50	208 × 175 × 34
631-1500	Blue	100	171 × 222 × 33
631-1501	Green	100	171 × 222 × 33
631-1502	Red	100	171 × 222 × 33
631-1503	Yellow	100	171 × 222 × 33
631-1504	White	100	171 × 222 × 33

Foam Lined

Catalogue Number	Colour	No. of Slides	Dimension [mm]
631-1522	Blue	25	141 × 88 × 35
631-1524	Red	25	141 × 88 × 35
631-1526	White	25	141 × 88 × 35
631-1517	Blue	100	171 × 222 × 33
631-1519	Red	100	171 × 222 × 33
631-1521	White	100	171 × 222 × 33

Quality System Compliance

Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements. inspected in according to AQL Level II 4.0 sampling plan per MIL-STD-1916.

Product Specifications

Material	Acrylonitrile-Butadiene Styrene.
Sterilization	Not applicable.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2021-07-08 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.

Storage Conditions Storage temperatures -20 °C to +60 °C.

Shelf Life 5 Years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.