

# Certificate of Conformity

<b>VWR Catalogue Number</b>	See List
<b>Description</b>	Zip-Seal Bags
<b>Country of Origin</b>	Manufactured in China
<b>Date of Issue (yyyy-mm-dd)</b>	2021-08-19

Single Pocket, Non-Sterile

Catalogue Number	Dimension [ mm ]
129-0293	60 × 40
129-0294	80 × 60
129-0309	100 × 70
129-0295	120 × 80
129-0296	150 × 100
129-0297	170 × 120
129-0298	220 × 160
129-0299	250 × 180
129-0306	300 × 200
129-0307	350 × 250
129-0308	400 × 300

Single Pocket with writing area, Non-Sterile

Catalogue Number	Dimension [ mm ]
129-0310	60 × 40
129-0311	80 × 60
129-0320	100 × 70
129-0312	120 × 80
129-0313	150 × 100
129-0314	170 × 120
129-0315	220 × 160

Catalogue Number	Dimension [ mm ]
129-0316	250 × 180
129-0317	300 × 200
129-0318	350 × 250
129-0319	400 × 300

Single Pocket with writing area, Sterile

Catalogue Number	Dimension [ mm ]
129-0325	300 × 200

Double Pocket, Non-Sterile

Catalogue Number	Dimension [ mm ]
129-0321	170 × 120
129-0322	220 × 160
129-0323	250 × 180
129-0324	300 × 200

### Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

### QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements. Inspected according to AQL level II 4.0 sampling plan per MIL-STD-105E.

Weld strength	:	≥ 15N/25 mm.
Tensile Strength	:	(Vertical >10MP/25 mm / Horizontal > 20MP/25 mm.
Specific Elongation	:	(Vertical ≥ 450% / Horizontal ≥ 250%.

## Product Specifications

<b>Material</b>	<p>Low Density Polyethylene (LDPE) C150Y + Linear Low-Density Polyethylene (LDPE) 218W (virgin).</p> <p>Monomers and additives used to produce the products are listed in List of Authorised Substances of Regulation EU/10/2011.</p> <p>The material grades conform to EU Directive 1935/2004 on materials and articles intended to come into contact with food.</p>
<b>Sterilization</b>	<p>Product labeled as sterile is gamma irradiated and dose released upon ISO11137 recommended practices in effect at the time of validation and meets a minimum requirement of 10<sup>-6</sup> SAL Which is not applicable for Non-Sterile products.</p>
<b>ATP Assay</b>	Not applicable.
<b>DNase &amp; RNase Free</b>	Not applicable.
<b>BSE/TSE</b>	<p>No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).</p>
<b>Cytotoxicity</b>	Not applicable.
<b>Non-Pyrogenic Statement</b>	Not applicable.
<b>Latex Statement</b>	This product is not made from natural rubber latex.
<b>BPA Statement</b>	Not applicable.
<b>DEHP Statement</b>	Not applicable.
<b>RoHS</b>	Not applicable.
<b>REACH Statement</b>	Not applicable.
<b>Storage Conditions</b>	Room temperature.
<b>Shelf Life</b>	Products labelled as sterile has shelf life period of 3 years which is not applicable for non-sterile products.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.  
This document has been produced electronically and is valid without a signature.