



EU DECLARATION OF CONFORMITY

Name and address of the manufacturer: / **Ningbo Transo-Medical Co., Ltd.**
Room 1608, No.18 Diamond Commercial plaza, No.409 Jiang'an Road, Jiangbei District, Ningbo City, Zhejiang Province, P.R. China
SRN: CN-MF-000009687

Registered trade name or mark:/ **None**

EU Authorized Representative:/ **Caretechion GmbH**
Niederrhein Straße, 71, Düsseldorf, Germany
info@ caretechion.com
SRN: DE-AR-000005946

We declare under our sole responsibility that

Name of the medical device: / **Cotton Swabs**

Product code:/ **Wooden swabs: TM1301**
Plastic swabs: TM1302

UMDN Code: 15066 Swabs, Cotton
EMDN Code: A1101, DEVICES FOR ADMINISTRATION, WITHDRAWAL AND COLLECTION, SAMPLE COLLECTION NEUTRAL SWABS

Intended purpose: / **The cotton swab is a tool for applying drug or disinfectant to disinfect the skin. The device does not contain drugs or disinfectants, for single-use only.**

Basic UDI-DI:/ **697461841TM1301JB (Wooden swabs)**
697461841TM1302JD (Plastic swabs)

of class: / **Rule 1, Class I**

According to annex VIII of Regulation (EU) 2017/745 /

Conformity assessment procedure:/ **Declare the conformity of the abovementioned products by issuing this EU Declaration of Conformity after drawing up the technical documentation set out in Annexes II and III of Regulation (EU) 2017/745 / according to Article 52(7) of Regulation (EU) 2017/745 /**

CS reference: / **ISO 15223-1:2021; ISO 14971:2019; ISO 10993-1:2018; ISO 10993-5:2009; ISO 10993-10:2013; ISO 11737-1:2018**

Meets the provisions of the Regulation EU 2017/745(MDR) which apply to it. The declaration is valid in connection with the "final inspection report" of the device. /

Ningbo China Nov 26th, 2021

Chenyun Ye, QM Department

Ort, Datum / Place, date /
Lieu, date / Luogo, data

Name und Funktion / Name and function /
Nom et fonction / Nome e funzione

Ningbo Transo-Medical Co., Ltd.
Cotton Swabs Finished Product Inspection Report

QS-FM-076 (C/1)

Client	Heinz Herenz Medizinalbedarf f GmbH	Art No.	1030128	Contract No.	22TM187 (62200913)	Lot No.	220901a
Description	wooden applicators 150x2.2mm small head			Quantity	2,680,000pcs	Date	20220916

Packing	Item	Inspection standards and requirements	Comformity	Defect description	Failed	Total failed	Qualified
	Size	According to TDS	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No
	Material	According to TDS	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0		
	Printing	Clear, visible and consistent with the artwork	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0		

Visual inspection	Item	Inspection standards and requirements	Sample size	Test Result			Qualified
				Defect description	Failed	Total failed	
	Foreign body	No foreign body	G2 AQL 2.5 1250 pcs Ac:21 Re:22	/	/	6	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No
	Surface rough	Free of damage and prick		wooden Stick damage	4		
	Shape	Head in good shape		/	/		
	Dirt	No yellowing and black spots		yellow	2		
	Other	No other visible defects		/	/		

Product dimension	Item	Inspection standards and requirements	Sample size	Test Result											Qualified	
				Average	Failed	Total failed										
	Length 1	According to TDS, 149~152mm	S2 AQL 4.0 13 pcs Ac: 1 Re: 2	150.5mm	0	0									<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	
	Length 2	According to TDS, 10~14mm		12.3mm	0											
	Dia. 1	According to TDS, 2.1~2.3mm		2.2mm	0											
	Dia. 2	According to TDS, 4~5mm		4.65mm	0											
	Data															
	Item	1	2	3	4	5	6	7	8	9	10	11	12	13	Test Instrument	
	Length 1	150	150	151	150	150	150	151	150	152	151	150	150	151	Ruler	
	Length 2	12	12	12	13	13	12	14	11	12	12	12	12	13	Ruler	
Dia. 1	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	Vernier caliper		
Dia. 2	4.6	4.6	4.6	4.8	4.8	4.8	4.8	4.6	4.4	4.4	4.6	4.8	4.6	Vernier caliper		

Performace	Item	Inspection standards and requirements	Sample size	Test Result			Qualified
				Defect description	Failed	Total failed	
	Head firmness	The cotton head will not fall off under a certain tension	S2 AQL 4.0 13 pcs Ac: 1 Re: 2	/	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No

Test Report	<input checked="" type="checkbox"/> Pass <input type="checkbox"/> Partial rework <input type="checkbox"/> Total rework <input type="checkbox"/> Recheck <input type="checkbox"/> Redo <input type="checkbox"/> Other_____						
Inspector	watson			Date		2022/9/16	

Comments:

Ningbo Transo-Medical Co., Ltd.
Cotton Swabs Finished Product Inspection Report

QS-FM-076 (C/1)

Client	Heinz Herenz Medizinalbedarf f GmbH	Art No.	1030818	Contract No.	22TM187 (62200913)	Lot No.	220901b
Description	Wooden applicator 200x3mm big head			Quantity	500,000pcs	Date	20220916

	Item	Inspection standards and requirements	Comformity	Defect description	Failed	Total failed	Qualified
Packing	Size	According to TDS	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No
	Material	According to TDS	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0		
	Printing	Clear, visible and consistent with the artwork	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0		

	Item	Inspection standards and requirements	Sample size	Test Result			Qualified
				Defect description	Failed	Total failed	
Visual inspection	Foreign body	No foreign body	G2 AQL 2.5 800 pcs Ac:21 Re:22	Dirty	2	4	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No
	Surface rough	Free of damage and prick		/	/		
	Shape	Head in good shape		/	/		
	Dirt	No yellowing and black spots		yellow/ black spots	2		
	Other	No other visible defects		/	/		

	Item	Inspection standards and requirements	Sample size	Test Result			Qualified									
				Average	Failed	Total failed										
Product dimension	Length 1	According to TDS,203~208mm	S2 AQL 4.0 13 pcs Ac: 1 Re: 2	205.1mm	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No									
	Length 2	According to TDS,27~32mm		29.4mm	0											
	Dia. 1	According to TDS,2.9~3.1mm		2.9mm	0											
	Dia. 2	According to TDS,10~11mm		10.5mm	0											
	Data															
		Item	1	2	3	4	5	6	7	8	9	10	11	12	13	Test Instrument
		Length 1	205	205	204	206	206	206	205	205	206	204	204	205	205	Ruler
		Length 2	30	28	30	29	29	29	28	28	29	30	30	30	32	Ruler
		Dia. 1	2.9	2.9	2.9	2.9	2.9	2.9	2.9	2.9	2.9	2.9	2.9	2.9	2.9	Vernier caliper
		Dia. 2	10.4	10.6	10.6	1.0.6	10.4	10.4	10.6	10.4	10.4	10.6	10.2	10.6	10.6	Vernier caliper

	Item	Inspection standards and requirements	Sample size	Test Result			Qualified
				Defect description	Failed	Total failed	
Performance	Head firmness	The cotton head will not fall off under a certain tension	S2 AQL 4.0 13 pcs Ac: 1 Re: 2	/	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No

Test Report	<input checked="" type="checkbox"/> Pass <input type="checkbox"/> Partial rework <input type="checkbox"/> Total rework <input type="checkbox"/> Recheck <input type="checkbox"/> Redo <input type="checkbox"/> Other_____
Inspector	Watson
	Date
	2022/9/16

Comments:

Ningbo Transo-Medical Co., Ltd.
Cotton Swabs Finished Product Inspection Report

QS-FM-076 (C/1)

Client	Heinz Herenz Medizinalbedarf GmbH	Art No.	1030918	Contract No.	22TM187 (62200913)	Lot No.	220901c
Description	Wooden applicator 250x3mm big head			Quantity	20,000pcs	Date	20220916

Packing	Item	Inspection standards and requirements	Comformity	Defect description	Failed	Total failed	Qualified
	Size	According to TDS	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No
	Material	According to TDS	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0		
	Printing	Clear, visible and consistent with the artwork	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0		

Visual inspection	Item	Inspection standards and requirements	Sample size	Test Result			Qualified
				Defect description	Failed	Total failed	
	Foreign body	No foreign body	G2 AQL 2.5 315 pcs Ac:14 Re:15	Dirty	2	4	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No
	Surface rough	Free of damage and prick		/	/		
	Shape	Head in good shape		/	/		
	Dirt	No yellowing and black spots		yellow/ black spots	2		
	Other	No other visible defects		/	/		

Product dimension	Item	Inspection standards and requirements	Sample size	Test Result											Qualified	
				Average	Failed	Total failed										
	Length 1	According to TDS,253~258mm	S2 AQL 4.0 8 pcs Ac: 1Re: 2	254.8mm	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No									
	Length 2	According to TDS,27~32mm		29.88mm	0											
	Dia. 1	According to TDS,2.9~3.1mm		2.9mm	0											
	Dia. 2	According to TDS,10~11mm		10.58mm	0											
	Data															
	Item	1	2	3	4	5	6	7	8	9	10	11	12	13	Test Instrument	
	Length 1	254	255	255	255	254	256	255	254							Ruler
	Length 2	31	29	28	30	30	30	30	31							Ruler
Dia. 1	2.9	2.9	2.9	2.9	2.9	2.9	2.9	2.9							Vernier caliper	
Dia. 2	10.4	10.6	10.6	10.6	10.4	10.8	10.8	10.4							Vernier caliper	

Performace	Item	Inspection standards and requirements	Sample size	Test Result			Qualified
				Defect description	Failed	Total failed	
	Head firmness	The cotton head will not fall off under a certain tension	S2 AQL 4.0 8pcs Ac: 1Re: 2	/	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No

Test Report	<input checked="" type="checkbox"/> Pass <input type="checkbox"/> Partial rework <input type="checkbox"/> Total rework <input type="checkbox"/> Recheck <input type="checkbox"/> Redo <input type="checkbox"/> Other_____
Inspector	<div style="display: flex; justify-content: space-between;"> Watson Date 2022/9/16 </div>

Comments:

Ningbo Transo-Medical Co., Ltd.
Cotton Swabs Finished Product Inspection Report

QS-FM-076 (C/1)

Client	Heinz Herenz Medizinalbedarf f GmbH	Art No.	1031018	Contract No.	22TM187 (62200913)	Lot No.	220901d
Description	Wooden applicator 300x3mm big head			Quantity	40,000pcs	Date	20220916

	Item	Inspection standards and requirements	Comformity	Defect description	Failed	Total failed	Qualified
Packing	Size	According to TDS	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No
	Material	According to TDS	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0		
	Printing	Clear, visible and consistent with the artwork	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0		

	Item	Inspection standards and requirements	Sample size	Test Result			Qualified
				Defect description	Failed	Total failed	
	Visual inspection	Foreign body	No foreign body	G2 AQL 2.5 500 pcs Ac:21 Re:22	/	/	2
Surface rough		Free of damage and prick	/		/		
Shape		Head in good shape	/		/		
Dirt		No yellowing and black spots	yellow/ black spots		2		
Other		No other visible defects	/		/		

	Item	Inspection standards and requirements	Sample size	Test Result			Qualified									
				Average	Failed	Total failed										
	Product dimension	Length 1	According to TDS,303~308mm	S2 AQL 4.0 13pcs Ac: 1Re: 2	305mm	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No								
Length 2		According to TDS,27~32mm	29.4mm		0											
Dia. 1		According to TDS,2.9~3.1mm	2.9mm		0											
Dia. 2		According to TDS,10~11mm	10.6mm		0											
Data																
		Item	1	2	3	4	5	6	7	8	9	10	11	12	13	Test Instrument
		Length 1	305	304	305	306	305	305	306	305	305	306	304	305	305	Ruler
		Length 2	28	28	29	29	28	30	31	30	30	29	29	30	31	Ruler
		Dia. 1	2.9	2.9	2.9	2.9	2.9	2.9	2.9	2.9	2.9	2.9	2.9	2.9	2.9	Vernier caliper
	Dia. 2	10.6	10.6	10.8	10.6	10.6	10.8	10.8	10.6	10.6	10.6	10.4	10.4	10.6	Vernier caliper	

	Item	Inspection standards and requirements	Sample size	Test Result			Qualified
				Defect description	Failed	Total failed	
Performace	Head firmness	The cotton head will not fall off under a certain tension	S2 AQL 4.0 13pcs Ac: 1Re: 2	/	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No

Test Report	<input checked="" type="checkbox"/> Pass <input type="checkbox"/> Partial rework <input type="checkbox"/> Total rework <input type="checkbox"/> Recheck <input type="checkbox"/> Redo <input type="checkbox"/> Other_____						
Inspector	Watson			Date		2022/9/16	
Comments:							

Ningbo Transo-Medical Co., Ltd.
Cotton Swabs Finished Product Inspection Report

QS-FM-076 (C/1)

Client	Heinz Herenz Medizinalbedar f GmbH	Art No.	1031118	Contract No.	22TM187 (62200913)	Lot No.	220901e
Description	Wooden applicator 150x2.2mm 9mm head			Quantity	20,000pcs	Date	20220916

	Item	Inspection standards and requirements	Comformity	Defect description	Failed	Total failed	Qualified
Packing	Size	According to TDS	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No
	Material	According to TDS	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0		
	Printing	Clear, visible and consistent with the artwork	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0		

Visual inspection	Item	Inspection standards and requirements	Sample size	Test Result			Qualified
				Defect description	Failed	Total failed	
	Foreign body	No foreign body	G2 AQL 2.5 315 pcs Ac:14 Re:15	/	/	5	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No
	Surface rough	Free of damage and prick		wooden Stick damage	3		
	Shape	Head in good shape		/	/		
	Dirt	No yellowing and black spots		yellow	2		
	Other	No other visible defects		/	/		

Product dimension	Item	Inspection standards and requirements	Sample size	Test Result			Qualified									
				Average	Failed	Total failed										
	Length 1	According to TDS,150~154mm	S2 AQL 4.0 8 pcs Ac: 1Re: 2	152mm	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No									
	Length 2	According to TDS,25~30mm		26.88mm	0											
	Dia. 1	According to TDS,2.1~2.3mm		2.2mm	0											
	Dia. 2	According to TDS,8~10mm		9.23mm	0											
	Data															
	Item	1	2	3	4	5	6	7	8	9	10	11	12	13	Test Instrument	
	Length 1	153	152	152	152	152	153	151	151						Ruler	
	Length 2	27	26	27	27	26	27	27	28						Ruler	
	Dia. 1	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2						Vernier caliper	
Dia. 2	9.6	9.4	9.4	9	9.4	9	9.2	8.8						Vernier caliper		

Perform	Item	Inspection standards and requirements	Sample size	Test Result			Qualified
				Defect description	Failed	Total failed	
	Head firmness	The cotton head will not fall off under a certain tension	S2 AQL 4.0 8 pcs Ac: 1Re: 2	/	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No

Test Report	<input checked="" type="checkbox"/> Pass <input type="checkbox"/> Partial rework <input type="checkbox"/> Total rework <input type="checkbox"/> Recheck <input type="checkbox"/> Redo <input type="checkbox"/> Other_____
Inspector	Watson
	Date
	2022/9/16

Comments:

Ningbo Transo-Medical Co., Ltd.
Cotton Swabs Finished Product Inspection Report

QS-FM-076 (C/1)

Client	Heinz Herenz Medizinalbedarf f GmbH	Art No.	1031128	Contract No.	22TM187 (62200913)	Lot No.	220901f
Description	plastic applicator 150x4mm big head			Quantity	100,000pcs	Date	20220916

Packing	Item	Inspection standards and requirements	Comformity	Defect description	Failed	Total failed	Qualified
	Size	According to TDS	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No
	Material	According to TDS	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0		
	Printing	Clear, visible and consistent with the artwork	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0		

Visual inspection	Item	Inspection standards and requirements	Sample size	Test Result			Qualified
				Defect description	Failed	Total failed	
	Foreign body	No foreign body	G2 AQL 2.5 500 pcs Ac:21 Re:22	/	/	2	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No
	Surface rough	Free of damage and prick		/	/		
	Shape	Head in good shape		/	/		
	Dirt	No yellowing and black spots		yellow/ black spots	2		
	Other	No other visible defects		/	/		

Product dimension	Item	Inspection standards and requirements	Sample size	Test Result											Qualified	
				Average	Failed	Total failed										
	Length 1	According to TDS,147~152mm	S2 AQL 4.0 13pcs Ac: 1Re: 2	150mm	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No									
	Length 2	According to TDS,27~32mm		29.6mm	0											
	Dia. 1	According to TDS,3.7~3.9mm		3.7mm	0											
	Dia. 2	According to TDS,10~11mm		10.6mm	0											
	Data															
	Item	1	2	3	4	5	6	7	8	9	10	11	12	13	Test Instrument	
	Length 1	150	149	149	150	150	149	150	151	150	150	151	151	150	Ruler	
	Length 2	30	30	29	28	28	30	31	30	30	29	29	30	31	Ruler	
	Dia. 1	3.7	3.7	3.7	3.7	3.7	3.7	3.7	3.7	3.7	3.7	3.7	3.7	3.7	Vernier caliper	
	Dia. 2	10.6	10.6	10.8	10.6	10.6	10.8	10.8	10.6	10.6	10.6	10.4	10.4	10.6	Vernier caliper	

Performace	Item	Inspection standards and requirements	Sample size	Test Result			Qualified
				Defect description	Failed	Total failed	
	Head firmness	The cotton head will not fall off under a certain tension	S2 AQL 4.0 13pcs Ac: 1Re: 2	/	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No

Test Report	<input checked="" type="checkbox"/> Pass <input type="checkbox"/> Partial rework <input type="checkbox"/> Total rework <input type="checkbox"/> Recheck <input type="checkbox"/> Redo <input type="checkbox"/> Other_____						
Inspector	Watson			Date		2022/9/16	

Comments:

Ningbo Transo-Medical Co., Ltd.
Cotton Swabs Finished Product Inspection Report

QS-FM-076 (C/1)

Client	Heinz Herenz Medizinalbedarf f GmbH	Art No.	1031618	Contract No.	22TM187 (62200913)	Lot No.	220901g
Description	Wooden applicator 200x3mm long small head			Quantity	20,000pcs	Date	20220916

Packing	Item	Inspection standards and requirements	Comformity	Defect description	Failed	Total failed	Qualified
	Size	According to TDS	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No
	Material	According to TDS	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0		
	Printing	Clear, visible and consistent with the artwork	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0		

Visual inspection	Item	Inspection standards and requirements	Sample size	Test Result			Qualified
				Defect description	Failed	Total failed	
	Foreign body	No foreign body	G2 AQL 2.5 315 pcs Ac:14 Re:15	/	/	1	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No
	Surface rough	Free of damage and prick		/	/		
	Shape	Head in good shape		/	/		
	Dirt	No yellowing and black spots		yellow	1		
	Other	No other visible defects		/	/		

Product dimension	Item	Inspection standards and requirements	Sample size	Test Result											Qualified	
				Average	Failed	Total failed										
	Length 1	According to TDS,201~204mm	S2 AQL 4.0 8pcs Ac:1 Re:2	202.5mm	0	0									<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	
	Length 2	According to TDS,27~32mm		29.38mm	0											
	Dia. 1	According to TDS,2.9~3.1mm		2.9mm	0											
	Dia. 2	According to TDS,5~6mm		5.7mm	0											
	Data															
	Item	1	2	3	4	5	6	7	8	9	10	11	12	13	Test Instrument	
	Length 1	203	203	203	204	201	202	202	202						Ruler	
	Length 2	30	30	29	29	30	29	29	29						Ruler	
Dia. 1	2.9	2.9	2.9	2.9	2.9	2.9	2.9	2.9						Vernier caliper		
Dia. 2	5.4	5.8	6	6	5.8	5.6	5.6	5.6						Vernier caliper		

Perform ance	Item	Inspection standards and requirements	Sample size	Test Result			Qualified
				Defect description	Failed	Total failed	
	Head firmness	The cotton head will not fall off under a certain tension	S2 AQL 4.0 8pcs Ac:1 Re:2	/	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No

Test Report	<input checked="" type="checkbox"/> Pass <input type="checkbox"/> Partial rework <input type="checkbox"/> Total rework <input type="checkbox"/> Recheck <input type="checkbox"/> Redo <input type="checkbox"/> Other_____						
Inspector	Watson			Date		2022/9/16	
Comments:							

Ningbo Transo-Medical Co., Ltd.
Cotton Swabs Finished Product Inspection Report

QS-FM-076 (C/1)

Client	Heinz Herenz Medizinalbedarf f GmbH	Art No.	1031718	Contract No.	22TM187 (62200913)	Lot No.	220901h
Description	Wooden applicator 250x3mm long small head			Quantity	50,000pcs	Date	20220916

	Item	Inspection standards and requirements	Comformity	Defect description	Failed	Total failed	Qualified
Packing	Size	According to TDS	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No
	Material	According to TDS	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0		
	Printing	Clear, visible and consistent with the artwork	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0		

	Item	Inspection standards and requirements	Sample size	Test Result			Qualified
				Defect description	Failed	Total failed	
Visual inspection	Foreign body	No foreign body	G2 AQL 2.5 500pcs Ac:21 Re:22	/	/	2	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No
	Surface rough	Free of damage and prick		/	/		
	Shape	Head in good shape		/	/		
	Dirt	No yellowing and black spots		yellow	2		
	Other	No other visible defects		/	/		

	Item	Inspection standards and requirements	Sample size	Test Result			Qualified								
				Average	Failed	Total failed									
Product dimension	Length 1	According to TDS,251~254mm	S2 AQL 4.0 13pcs Ac:1 Re:2	251.9mm	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No								
	Length 2	According to TDS,27~32mm		28.9mm	0										
	Dia. 1	According to TDS,2.9~3.1mm		3mm	0										
	Dia. 2	According to TDS,5~6mm		5.9mm	0										
	Data														
	Item	1	2	3	4	5	6	7	8	9	10	11	12	13	Test Instrument
	Length 1	253	254	252	252	251	251	251	252	252	252	251	252	252	Ruler
	Length 2	28	29	29	29	30	28	29	29	29	29	29	30	28	Ruler
	Dia. 1	3	3	3	3	3	3	3	3	3	3	3	3	3	Vernier caliper
Dia. 2	5.4	5.8	6	6	5.8	5.6	5.6	5.6	6.2	6.2	6.4	6.2	6	Vernier caliper	

	Item	Inspection standards and requirements	Sample size	Test Result			Qualified
				Defect description	Failed	Total failed	
Perform ance	Head firmness	The cotton head will not fall off under a certain tension	S2 AQL 4.0 13pcs Ac:1 Re:2	/	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No

Test Report	<input checked="" type="checkbox"/> Pass <input type="checkbox"/> Partial rework <input type="checkbox"/> Total rework <input type="checkbox"/> Recheck <input type="checkbox"/> Redo <input type="checkbox"/> Other_____						
Inspector	Watson			Date		2022/9/16	
Comments:							

Ningbo Transo-Medical Co., Ltd.

QS-FM-076 (C/1)

Cotton Swabs Finished Product Inspection Report

Client		Heinz Herenz Medizinalbedarf f GmbH	Art No.	1031918	Contract No.	22TM187 (62200913)	Lot No.	220901q							
Description		wooden applicators 230x2.2mm small head			Quantity	1,000,000pcs	Date	20220916							
Packing	Item	Inspection standards and requirements	Comformity	Defect description	Failed	Total failed	Qualified								
	Size	According to TDS	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No								
	Material	According to TDS	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0										
	Printing	Clear, visible and consistent with the artwork	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0										
Visual inspection	Item	Inspection standards and requirements	Sample size	Test Result			Qualified								
				Defect description	Failed	Total failed									
	Foreign body	No foreign body	G2 AQL 2.5 1250 pcs Ac:21 Re:22	/	/	4	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No								
	Surface rough	Free of damage and prick		wooden Stick damage	3										
	Shape	Head in good shape		/	/										
	Dirt	No yellowing and black spots		yellow	1										
	Other	No other visible defects		/	/										
Product dimension	Item	Inspection standards and requirements	Sample size	Test Result			Qualified								
				Average	Failed	Total failed									
	Length 1	According to TDS,229~232mm	S2 AQL 4.0 13 pcs Ac: 1 Re: 2	230.1mm	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No								
	Length 2	According to TDS,10~14mm		12.3mm	0										
	Dia. 1	According to TDS,2.1~2.3mm		2.2mm	0										
	Dia. 2	According to TDS,4~5mm		4.6mmmm	0										
	Data														
	Item	1	2	3	4	5	6	7	8	9	10	11	12	13	Test Instrument
	Length 1	230	230	229	230	230	230	229	230	230	231	231	231	230	Ruler
	Length 2	12	12	13	12	12	12	12	13	13	13	12	12	12	Ruler
	Dia. 1	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	Vernier caliper
Dia. 2	4.6	4.4	4.6	4.4	4.6	4.6	4.6	4.6	4.6	4.8	4.8	4.6	4.4	Vernier caliper	
Perform ance	Item	Inspection standards and requirements	Sample size	Test Result			Qualified								
				Defect description	Failed	Total failed									
	Head firmness	The cotton head will not fall off under a certain tension	S2 AQL 4.0 13 pcs Ac: 1 Re: 2	/	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No								
Test Report		<input checked="" type="checkbox"/> Pass <input type="checkbox"/> Partial rework <input type="checkbox"/> Total rework <input type="checkbox"/> Recheck <input type="checkbox"/> Redo <input type="checkbox"/> Other_____													
Inspector		Watson			Date		2022/9/16								
Comments:															

Ningbo Transo-Medical Co., Ltd.
Cotton Swabs Finished Product Inspection Report

QS-FM-076 (C/1)

Client	Heinz Herenz Medizinalbedarf f GmbH	Art No.	1032018	Contract No.	22TM187 (62200913)	Lot No.	220901i
Description	Wooden applicator 75x2.2mm double small head			Quantity	500,000pcs	Date	20220916

Packing	Item	Inspection standards and requirements	Comformity	Defect description	Failed	Total failed	Qualified
	Size	According to TDS	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No
	Material	According to TDS	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0		
	Printing	Clear, visible and consistent with the artwork	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0		

Visual inspection	Item	Inspection standards and requirements	Sample size	Test Result			Qualified
				Defect description	Failed	Total failed	
	Foreign body	No foreign body	G2 AQL 2.5 800 pcs Ac:21 Re:22	/	/	2	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No
	Surface rough	Free of damage and prick		/	/		
	Shape	Head in good shape		/	/		
	Dirt	No yellowing and black spots		yellow	2		
	Other	No other visible defects		/	/		

Product dimension	Item	Inspection standards and requirements	Sample size	Test Result			Qualified								
				Average	Failed	Total failed									
	Length 1	According to TDS, 75~78mm	S2 AQL 4.0 13 pcs Ac: 1 Re: 2	76.3mm	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No								
	Length 2	According to TDS, 10~14mm		12.7mm	0										
	Dia. 1	According to TDS, 2.1~2.3mm		2.2mm	0										
	Dia. 2	According to TDS, 4~5mm		4.5mm	0										
	Data														
	Item	1	2	3	4	5	6	7	8	9	10	11	12	13	Test Instrument
	Length 1	77	75	75	77	77	78	77	75	77	76	76	76	76	Ruler
	Length 2	13	12	13	13	14	12	13	12	13	14	12	12	12	Ruler
Dia. 1	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	Vernier caliper	
Dia. 2	4.6	4.2	4.6	4.4	4.6	4.6	4.4	4.4	4.4	4.6	4.5	4.6	4.6	Vernier caliper	

Perform ance	Item	Inspection standards and requirements	Sample size	Test Result			Qualified
				Defect description	Failed	Total failed	
	Head firmness	The cotton head will not fall off under a certain tension	S2 AQL 4.0 13 pcs Ac: 1 Re: 2	/	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No

Test Report	<input checked="" type="checkbox"/> Pass <input type="checkbox"/> Partial rework <input type="checkbox"/> Total rework <input type="checkbox"/> Recheck <input type="checkbox"/> Redo <input type="checkbox"/> Other_____
Inspector	Watson
Date	2022/9/16

Comments:

Ningbo Transo-Medical Co., Ltd.
Cotton Swabs Finished Product Inspection Report

QS-FM-076 (C/1)

Client	Heinz Herenz Medizinalbedarf f GmbH	Art No.	1033001	Contract No.	22TM187 (62200913)	Lot No.	220901p
Description	wooden applicators 200x2.2mm small head			Quantity	60,000pcs	Date	20220916

	Item	Inspection standards and requirements	Comformity	Defect description	Failed	Total failed	Qualified
Packing	Size	According to TDS	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No
	Material	According to TDS	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0		
	Printing	Clear, visible and consistent with the artwork	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	/	0		

Visual inspection	Item	Inspection standards and requirements	Sample size	Test Result			Qualified
				Defect description	Failed	Total failed	
	Foreign body	No foreign body	G2 AQL 2.5 500 pcs Ac:21 Re:22	/	/	2	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No
	Surface rough	Free of damage and prick		/	/		
	Shape	Head in good shape		/	/		
	Dirt	No yellowing and black spots		yellow	2		
	Other	No other visible defects		/	/		

Product dimension	Item	Inspection standards and requirements	Sample size	Test Result											Qualified
				Average	Failed	Total failed									
	Length 1	According to TDS, 201~204mm	S2 AQL 4.0 13 pcs Ac: 1 Re: 2	202.5mm	0	0									<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No
	Length 2	According to TDS, 10~14mm		12.6mm	0										
	Dia. 1	According to TDS, 2.1~2.3mm		2.2mm	0										
	Dia. 2	According to TDS, 4~5mm		4.5mm	0										
	Data														
	Item	1	2	3	4	5	6	7	8	9	10	11	12	13	Test Instrument
	Length 1	203	203	201	202	202	202	203	202	203	203	203	203	203	Ruler
	Length 2	14	13	12	12	13	12	12	12	12	13	13	13	13	Ruler
	Dia. 1	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	Vernier caliper
	Dia. 2	4.4	4.6	4.6	4.4	4.4	4.6	4.2	4.6	4.4	4.6	4.6	4.4	4.4	Vernier caliper

Perform	Item	Inspection standards and requirements	Sample size	Test Result			Qualified
				Defect description	Failed	Total failed	
	Head firmness	The cotton head will not fall off under a certain tension	S2 AQL 4.0 13 pcs Ac: 1 Re: 2	/	0	0	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No

Test Report	<input checked="" type="checkbox"/> Pass <input type="checkbox"/> Partial rework <input type="checkbox"/> Total rework <input type="checkbox"/> Recheck <input type="checkbox"/> Redo <input type="checkbox"/> Other_____						
Inspector	Watson			Date		2022/9/16	

Comments: