

Gebrauchsfertige Lösung		1.000 ml = 1 l	2.000 ml = 2 l	3.000 ml = 3 l	4.000 ml = 4 l	5.000 ml = 5 l	8.000 ml = 8 l	10.000 ml = 10 l
0,5 % ige Lösung	Wasser (ml)	995,0	1.990,0	2.985,0	3.980,0	4.975,0	7.960,0	9.950,0
	Desinfektionslösung (ml)	5,0	10,0	15,0	20,0	25,0	40,0	50,0
1,0 % ige Lösung	Wasser (ml)	990,0	1.980,0	2.970,0	3.960,0	4.950,0	7.920,0	9.900,0
	Desinfektionslösung (ml)	10,0	20,0	30,0	40,0	50,0	80,0	100,0
1,5 % ige Lösung	Wasser (ml)	985,0	1.970,0	2.955,0	3.940,0	4.925,0	7.880,0	9.850,0
	Desinfektionslösung (ml)	15,0	30,0	45,0	60,0	75,0	120,0	150,0
2,0 % ige Lösung	Wasser (ml)	980,0	1.960,0	2.940,0	3.920,0	4.900,0	7.840,0	9.800,0
	Desinfektionslösung (ml)	20,0	40,0	60,0	80,0	100,0	160,0	200,0
3,0 % ige Lösung	Wasser (ml)	970,0	1.940,0	2.910,0	3.880,0	4.850,0	7.760,0	9.700,0
	Desinfektionslösung (ml)	30,0	60,0	90,0	120,0	150,0	240,0	300,0
4,0 % ige Lösung	Wasser (ml)	960,0	1.920,0	2.880,0	3.840,0	4.800,0	7.680,0	9.600,0
	Desinfektionslösung (ml)	40,0	80,0	120,0	160,0	200,0	320,0	400,0
5,0 % ige Lösung	Wasser (ml)	950,0	1.900,0	2.850,0	3.800,0	4.750,0	7.600,0	9.500,0
	Desinfektionslösung (ml)	50,0	100,0	150,0	200,0	250,0	400,0	500,0
7,5 % ige Lösung	Wasser (ml)	925,0	1.850,0	2.775,0	3.700,0	4.625,0	7.400,0	9.250,0
	Desinfektionslösung (ml)	75,0	150,0	225,0	300,0	375,0	600,0	750,0